



microSystem[®] - Dispenser MKII

The microSystem Dispenser with 25ml cartridges offers significant benefits over the 50ml systems currently available.

31 % less waste – saves money

28 % less weight – gentle and precise application

26 % smaller profile – does not restrict the intraoral view

No blocked cartridges.

No cross contamination.

The patented Coltène microSystem cartridges and mixing tips avoid contamination. Outflow apertures are set apart and are combined with a safety reservoir which avoids undesirable cross-contamination of material flowing back. A premature contact of both materials and subsequent undesirable polymerisation is prevented with the mixing tips.

31% less waste – saves money


microSystem minimises waste material with each impression with its specially developed mixing tips.



Gentle and precise application.

- The ergonomic microSystem Dispenser is perfectly suited to both small and large hands.
- The Dispenser is easy to clean and autoclavable for complete infection control.
- The accurate and smooth extrusion mechanism facilitates precise intra-oral application of the impression material

The picture below clearly shows the difference in size; the real tactile advantage will be immediately apparent during use.



The perfect combination with AFFINIS®

- The compact size of the Coltène® microSystem Dispenser provides relaxed, comfortable, impression taking.
- Using microSystem you can choose AFFINIS light and regular body. The “wash” impression materials are compatible with AFFINIS heavy body, MonoBody and putties.
- The unique surface affinity and the optimized characteristics of AFFINIS light- and regular body, in combination with the microSystem Dispenser MKII, are perfectly adjusted to each other.

AFFINIS® influences your success. All down the line.



Order Information

6500	AFFINIS light body microSystem 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (small), 20 pcs Oral Tips
6600	AFFINIS fast light body microSystem 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (small), 20 pcs Oral Tips
6510	AFFINIS regular body microSystem 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (medium), 20 pcs Oral Tips
6610	AFFINIS fast regular body microSystem 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (medium), 20 pcs Oral Tips

6779	AFFINIS®PRECIOUS light body microSystem 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tip Green, 20 pcs Oral Tips
6780	AFFINIS®PRECIOUS regular body microSystem 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tip Clear, 20 pcs Oral Tips

Accessories: microSystem®

6270	Dispenser microSystem	1 pcs
6214	Mixing Tip Clear	40 pcs
6215	Mixing Tip Clear	100 pcs
6202	Mixing Tip Green	40 pcs
6205	Mixing Tip Green	100 pcs
6210	Oral Tips	100 pcs

6530	AFFINIS putty soft Base and Catalyst, each 300 ml
6531	AFFINIS fast putty soft Base and Catalyst, each 300 ml
6535	AFFINIS putty super soft Base and Catalyst, each 300 ml
6520	AFFINIS heavy body 75 ml 2 × 75 ml Cartridges, 8 pcs Mixing Tips (turquoise/large)
6620	AFFINIS fast heavy body 75 ml 2 × 75 ml Cartridges, 8 pcs Mixing Tips (turquoise/large)
6547	AFFINIS MonoBody 75 ml 2 × 75 ml Cartridges, 8 pcs Mixing Tips (turquoise/large)

Accessories: 50/75 System

6900	Dispenser 50/75 ml	1 pc
6552	Mixing Tip Turquoise	40 pcs

6487	AFFINIS System 360 heavy body 2 × 300 ml Refill package Base, 2 × 62 ml Refill package Catalyst
6485	AFFINIS System 360 heavy body Starter Kit 1 × 300 ml Base, 1 × 62 ml Catalyst, 1 × Supporting Cartridges, 10 × Dynamic Mixing Tips
6482	AFFINIS MonoBody System 2 × 300 ml Refill package Base, 2 × 62 ml Refill package Catalyst
6480	AFFINIS System 360 MonoBody Starter Kit 1 × 300 ml Base, 1 × 62 ml Catalyst, 1 × Supporting Cartridges, 10 × Dynamic Mixing Tips

Accessories: System 360

6160	Supporting Cartridges	1 pc
6161	Dynamic Mixing Tip Green	40 pcs

www.coltenewhaledent.com
© 2010 Coltène/Whaledent

Coltène/Whaledent Ltd.
The President Suite-A
Kendal House, Burgess Hill
Victoria Way
West Sussex, RH15 9NF/U.K.
Tel. +44 (0)1444 235486
Fax +44 (0)1444 870640
info@coltenewhaledent.co.uk

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten/Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info@coltenewhaledent.ch

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223/USA
Tel. +1 330 916 8800
Fax +1 330 916 7077
info@coltenewhaledent.com

**coltène
whaledent**